

PRO-VET DUAL ROTOR CENTRIFUGE 2mL, 15mL & Haematocrit



Display indicative only.

Dims HWD235 x 235 x 350mm
Weight.....12.8 Kg (without rotor)
Power.....125 watts

Program	Speed	Time
Program 5401		
Haematocrit (Pcv)	12,000Rpm	300Secs
Serum 1	6,000Rpm	600Secs
Serum 2	5,000Rpm	300Secs
Urine	3,500Rpm	300Secs
Program 5408		
Serum 1	6,000Rpm	180Secs
Serum 2	6,000Rpm	240Secs
Serum 3	5,000Rpm	300Secs
Urine	3,500Rpm	300Secs

Serum (15ml max) & Haematocrit Rotor



Rotors

Rotor type	24 x Capillary & 12 x 2ml	8 x 15ml
Tube size max	2 x 75mm & 11 x 40mm	17 x 120mm
Minimum speed	500 Rpm	500 Rpm
Maximum speed	12,000 Rpm	6000 Rpm
Maximum Rcf (G)	13500	3500
Radius max cms	8.5	8.50
Sample tube angle (°)	0 & 60	30
Acceleration time (secs)	30	30
Deceleration time (secs)	30	30
Autoclaveable (frequency)	121°C (10)	121°C (20)

Accessories



PCV reader



Sealing clay