

VET-800 INFUSION PUMP



- Light weight compact design (only 1.8kg)
- Large 3.2" display and easy operation
- Adjustable occlusion pressure settings (3 levels)
- Adjustable ultrasonic Air In Line detection (3 levels)
- Supports a wide variety of 20 drop/ml infusion sets
- Optional fluid warmer

Flow rate range:	0.1 - 1200 ml/h
VTBI:	1 - 9999 ml
KVO rate:	0.1 - 10 ml/h
BOLUS:	Default 600 ml/h (adjustable)
Purge rate:	600 ml/h
Accuracy:	≤ ± 5% (after callibration)
Size:	167 x 125 x 115 mm
Weight:	1.8kg
Audio and visual alarms:	Air in line, Door open, Occlusion, Finished, etc
Battery:	4 hr operation at 125 ml/h after full charge