

SURGILITE

Ceiling Mounted Examination light



Features

- TGA registered, ARTG#130080
- Suitable for installation into cardiac protected electrical areas.
- Durable stainless steel construction.
- Filter to remove heat from light spot.
- No exposed glass components.
- Shatter proof polycarbonate front filter
- Continuous rotation at main pivot points to aid positioning of lights.
- Fully enclosed electrical wiring.
- Switch located on light head.
- Large range of vertical movement.
- External Safety Extra Low Voltage (SELV) Transformer.

Technical Specifications.

Light Output at 0.5m (max)	50,000Lux.
Light Output at 1.0m (max)	24,000Lux.
Size of Light Spot at 0.5m	10cm
Size of Light Spot at 1.0m	12cm
Colour Temperature (Kelvin)	3400 K
Total Power Consumption	45VA
Bulb Type	OSRAM 12V 35W 4° P/N 41832 SSP
Transformer Type	ATCO TM50-2
Range of Vertical Movement	120cm
Standard Ceiling Heights (others available on request)	2500, 2700, 3000, 3300mm
Warranty Period	3 Years

