


NUWEIGH

VETERINARY SCALES

CHR392

NUWEIGH'S Veterinary Scales with large stainless steel platform uses four precision loadcells for high accuracy readings.



Large 26mm backlit display features a magnetic wall mount bracket to attach indicator remotely from weighing platform




FUNCTIONS



FEATURES

- Four precision loadcells
- Auto OFF function for power saving
- Animal hold function for effortless recordings
- Adjustable feet
- Reads in both kgs, oz and lbs
- 9V power adaptor or 4 x AA battery

DIMENSIONS



Capacity	Divisions	Part No.
250kg	100g	#36617

