

NOVO B+ 23L B-CLASS AUTOCLAVE

- Cycles for porous, packed and unpacked loads
- Suitable for solid and hollow instruments
- Supplied with circular rack system, trays, and manual
- Chamber dimensions: 245mm x 470mm Deep
- Outer dimensions: 670 x 450 x 400mm
- Sterilizing temperature: 121°C and 134°C

Program List

Program	Temperature	Pressure	Pre-Vacuum	SterilizationTime	DryingTime
UNPACKED	134°C	210kPa	1	4 min	9 min
PACKED	134°C	210kPa	3	5 min	9 min
PRION	134°C	210kPa	3	18 min	9 min
POROUS	121°C	110kPa	3	20 min	18 min

